

Keep Going<sup>sm</sup> with



The OxyGo<sup>®</sup> family of Portable Oxygen Concentrators gives patients reliable oxygen—without bulky tanks!

- Airline approved
- Up to 6 intelligent pulse flow settings
- Lightweight, quiet and discreet
- Easy-to-replace batteries lasting up to 13 hours
- Now with the **NEW “My OxyGo” App!**\*



[www.oxygo.life](http://www.oxygo.life)

\* Bluetooth available on OxyGo FIT POCs produced after 2/19 and OxyGo NEXT.

# Is OxyGo® right for your patients?

**ABSOLUTELY!** Thousands of patients love the OxyGo family of Portable Oxygen Concentrators, and so will yours. With their compact size, light weight, and long battery life, you can clearly see how OxyGo outperforms all other concentrators on the market. What makes OxyGo POCs so unique is their extremely high reliability rate, outstanding warranty, and innovative design that reduces unit down-time.



## OxyGo NEXT™

7.19" x 3.26" x 8.15" 4.7 lbs.

**OxyGo NEXT is the next generation of portable oxygen.** Its ground-breaking size and battery life have never been seen before at this high of an oxygen dose, and it is setting the standard for portable oxygen concentrators. OxyGo NEXT has user-replaceable batteries, sieve beds and an easy to read LCD display with simple controls. And it can be charged at home or on the go—giving your patients the freedom to embrace what's NEXT!

# 1400-3000

# 1400-3000-16 (system with double cell battery)

- 6 Pulse Flow Settings (1260 ml/min)
- Up to 13 hour battery life\*
- Only 38 dB. at setting 2
- Included 5-year warranty on the unit; 1-year warranty on sieve beds & battery
- Bluetooth-enabled units connect to your patient's mobile device with the **My OxyGo App.\*\*** You and your patient can troubleshoot, check battery life, oxygen purity and more with the touch of a button!



## OxyGo FIT™

5.9" x 2.7" x 7.2" 2.8 lbs.

OxyGo FIT weighs just 2.8 pounds and has 3 settings, with 210 ml of oxygen per setting. The OxyGo FIT is well constructed, will operate for almost 3 hours using a standard, user-replaceable battery, and sieve beds. It also has an optional larger battery, which will power the unit for up to 5 hours\* and adds less than one-half pound in weight.

# 1400-2000

# 1400-2000-8 (system with double cell battery)

- 3 Pulse Flow Settings (630 ml/min)
- Up to 5 hour battery life\*
- Only 40 dB. at setting 2
- Included 3-year warranty on the unit; 1-year warranty on sieve beds & battery
- Bluetooth-enabled units connect to your patient's mobile device with the **My OxyGo App.\*\*** You and your patient can troubleshoot, check battery life, oxygen purity and more with the touch of a button!



### STAY CONNECTED TO YOUR PATIENTS

#### Provider Portal

- **Real Time Monitoring:** Concentrator health checks, Error troubleshooting, Software update pushes
- **Routine Maintenance:** Collect usage feedback to drive improvements, preventative maintenance
- **Geo-locating:** Locate assets, location-specific push notifications, promotional marketing

#### My OxyGo App

- Check battery life
- Monitor column life status
- Easy access to FAQs, troubleshooting, user manuals... and more!



Lightweight

Ultra-Light

	OxyGo NEXT™	CAIRE® Freestyle Comfort®	Respironics SimplyGo Mini®	Invacare® Platinum Mobile™	GCE Zen-O Lite	O <sub>2</sub> Concepts® Oxlife Freedom™	Drive POC® iGo 2™	OxyGo FIT™
Brand rated #1 by Industry Experts *	✓	✗	✗	✗	✗	✗	✗	✓
Brand Most Requested by Patients	✓	✗	✗	✗	✗	✗	✗	✓
Bluetooth Enabled	✓	✓	✗	✓	✗	✗	✗	✓
Rated IP22 Water Resistant	✓	✓	✓	✓	✓	✓	✓	✓
O <sub>2</sub> Capacity	1260 ml / min	1050 ml / min	1000 ml / min	880 ml / min	1050 ml / min	800 ml / min	1114 ml / min	630 ml / min
Can Be Used 24/7	✓	✗	✗	✗	✗	✗	✗	✓
User-Replaceable Batteries	✓	✓	✓	✓	✓	✓	✓	✓
Intelligent Delivery Technology	✓	✗	✗	✗	✗	✗	✗	✓
Flow Settings	1-6	1-5	1-5	1-4	1-5	1-5	1-5	1-3
Weight	4.7 lbs.	5 lbs.	5 lbs.	4.9 lbs.	5.5 lbs.	5.8 lbs.	5 lbs.	2.8 lbs.
Battery Duration	UP TO 6.5 hours single battery 13 hours double battery	UP TO 4.25 hours single battery	UP TO 4.5 hours	UP TO 5 hours	UP TO 4 hours	UP TO 4 hours	UP TO 3.5 hours	UP TO 2.7 hours single battery 5.0 hours double battery
Sound Level	38 dB.	43 dB.	42 dB.	40 dB.	44 dB.	56 dB.	41 dB.	40 dB.
Warranty	5 years unit, 1 year sieve beds & battery	3 years unit, 1 year battery & power supply	3 years unit, 1 year sieve beds	3 years unit, 2 years compressor, 1 yr. sieve beds	3 years unit or 15,000 hours	3 years unit, 1 year sieve beds & battery	5 years unit, 3 years compressor, 1 yr. sieve beds & battery	3 years unit, 1 year sieve beds & battery
Expected Compressor Life	20,000 hours continuous	Not Published	Not Published	4,380 hours continuous	Not Published	Not Published	Not Published	25,000 hours continuous

Information from manufacturer manuals as of July 2019.

\* from independent, third-party *Needham Report*

## OxyGo keeps patients going— plus local provider support to ensure they get the best care.

As a family-owned and operated company, we believe the highest quality of care comes from the highest quality product—plus close interaction and monitoring. We provide the highest quality product, but realize that local providers will be able to offer patients much more dedicated attention and service. That is why we do not sell OxyGo or its accessories directly to patients. This "teamwork" mentality clearly makes OxyGo the right choice to keep patients' oxygen needs satisfied, along with giving them the high standard of care they deserve. OxyGo Portable Oxygen Concentrators are the same patient-demanded, reliable concentrators as seen on TV—backed by the care and attention of the patient's local provider.



# Accessories & Replacement Parts

Description	OxyGo 5 Setting Part Number	OxyGo FIT Part Number	OxyGo NEXT Part Number
Single Lithium-ion rechargeable battery	1400-1010-8	1400-2010-4	1400-3010-8
Double Lithium-ion rechargeable battery	1400-1010-16	1400-2010-8	1400-3010-16
External / Desktop Battery Charger (includes charger, power supply and AC power cord)	1400-1030	1400-2030	1400-3030
Home / Office AC Power Supply (includes power supply and AC power cord)	1400-1040	1400-2040	1400-3040
Car / Boat / RV DC Power Supply (for use with OxyGo 5 Setting or OxyGo FIT)	1400-1050	1400-1050	1400-1050
Carrying Bag / Case	1170-1410	1170-2410	1170-3410
Carry Strap (only)	1170-1415	1170-2415	N/A
Backpack	1170-1420	1170-2420	1170-3420
Slim Backpack	N/A	1170-2440	1170-3441
Accessory Bag	1170-1445	1170-1445	1170-1445
Cart	1170-1480	N/A	N/A

## REPLACEMENT PARTS

Intake Particle Filters	1400-1310 & 1400-1311	1400-2310	1400-3310
Outlet Barb Spanner Wrench (qty. 1)	1400-1321	1400-1321	1400-1321
Oxygen Product Filter Replacement Kit (includes spanner wrench and 2 filters)	1400-1322	1400-2322	1400-2322
Molecular Sieve Beds (pair)	1400-1060 & 1400-1061	1400-2060	1400-3060
Fuse for DC Power Supply (#1400-1050)	1400-1051	1400-1051	1400-1051
User Manual	1400-MANUAL	1400-2421	1400-3421

## Prescribing OxyGo is EASY!

OxyGo is covered by Medicare and most insurances. To ensure your patients get an OxyGo POC, simply prescribe:

**“OxyGo NEXT or OxyGo FIT, no substitutions.”**

FOLLOW US:



[www.oxygo.life](http://www.oxygo.life)