



Forbes
HEALTH
REVIEWS

BUZZ
AROUND
Carry On[®]

AND NAMES IT AMONG
**THE BEST LIGHTWEIGHT
MOBILITY SCOOTERS OF 2022!**

 **GOLDEN**
YOUR LIFE IN MOTION



A MESSAGE FROM RICH GOLDEN, CEO

Golden is known for comfort, quality and innovation. Founded in 1985, we manufacture mobility products that push the boundaries of expectation, providing extremely comfortable, highly dependable scooters and power wheelchairs with superior designs to last for years to come.

Through our products, we help provide people who have mobility issues an avenue to pursue their own independence allowing them to live their best lives and find happiness. Our scooters and power wheelchairs offer excellent maneuverability and performance to get you where you want to go in comfort and style. Explore our line of mobility products and choose the one that's perfect for you! Plus, all models are backed by the best warranties in the industry!

We are proud to manufacture America's favorite scooters and power wheelchairs - for your family, and ours.

A handwritten signature in black ink that reads "Rich Golden".

OUR VISION...

To build the most innovative and comfortable products that put you and your family's lives in motion.

Golden scooters and power wheelchairs are Class II Medical Devices.



BUZZTM AROUND SERIES



Going Places

BUZZAROUND SERIES

Relax & Ride in Style!

Our scooters and power wheelchairs offer excellent maneuverability and performance to get you where you want to go in comfort and style. Explore our line of mobility products and choose the one that's perfect for you!

*Highly Dependable • Maximum Maneuverability
Unparalleled Stability • Greatest Functionality
Superior Designs • First Class Luxury*



**RATED BEST IN CLASS
RANGE & BEST IN CLASS
WEIGHT CAPACITY!**

Full-size comfort comes with the convenience of true portability on the Buzzaround LX, along with standard features that include a stylish captain's seat, an infinitely adjustable tiller at your fingertips, and a USB charging port right under the dashboard. With a 375 lb. weight capacity, comfort suspension springs, and up to 18 mile range on two U1 batteries, the Buzzaround LX offers ample foot room to take you where you want to go in comfort and style!

BUZZAROUND LX FEATURES:

- Stylish Captain's Seat
- 375 lb. Weight Capacity
- Infinitely Adjustable Tiller with EZ Reach Adjustment Lever
- Up to 18 Mile Range with two U1 Batteries
- Ample Foot Room
- USB Charging Port
- Comfort Spring Suspension for a Smoother Ride
- Red & Blue Color Shrouds Standard

OPTIONAL Y CABLE BATTERY CHARGING HARNESS!



Optional color shrouds available.
*Upcharge applies.



GB119



GB149

BUZZTM AROUND EX

EXTREME
FULL-SIZE PORTABLE

**RATED BEST IN
CLASS RANGE!**

The Golden BuzzAround Extreme Scooter features the portability and easy disassembly of a lightweight scooter with the amazing benefits of a full size scooter. Both the three and four-wheel BuzzAround Extreme Scooters come standard with 2-U1 batteries - industry leading and the largest battery in this category - for maximum operating range up to 18 miles. The Extreme also features a brand new ultrabright LED headlight and taillight, and comfort suspension springs for the smoothest ride. Go the extra mile on a BuzzAround Extreme Scooter!

BUZZAROUND EX FEATURES:

- 350 lb. weight capacity
- 4" ground clearance
- 2-U1 batteries for maximum range
- Ultrabright LED headlight and taillight
- Comfort suspension springs
- Operating range up to 18 miles per charge



OPTIONAL Y CABLE BATTERY CHARGING HARNESS!

The Y Cable Charging Harness conveniently connects directly from the electrical outlet to each battery box.



OPTIONAL BATTERY DOCKING STATION

The Battery Docking Station for models GB118, GB148, GB119 and GB149 is portable. It connects and charges the U1 batteries. Perfect for charging batteries off the scooter, especially when visiting family or friends, or at a vacation destination. Enjoy being on the go!



OPTIONAL BUMPER CONVERSION KIT

For easy vertical storage. Can be ordered as a KIT or pre-installed on models GB118, GB148, GB119 & GB149.



Optional color shrouds available.
*Upcharge applies.

BUZZAROUND SERIES



GB118



GB148

SCAN FOR VIDEO



BUZZXLSM
AROUND **XLSM**
HD

**Upgrade with NEW
XR Extended Range
22 AH Battery Option!***



GOLDEN



Jeannie D. | Google Review

"I absolutely love my Golden scooter because it comes apart into 5 pieces that fits easily into my little car's trunk. That gives me the freedom to be able to go EVERYWHERE! My scooter is a much-needed medical device that has thankfully given me back my life. I am eternally grateful that I have it."

The next generation of the popular Buzzaround scooter series is now available! The Buzzaround XL and heavy-duty XLSHD models offer improved drivability and standard features to make your riding experience better than ever. They feature a larger wraparound Delta tiller, higher ground clearance, and an LED lighting package. A new stylish, contoured stadium style seat offers more comfort and increased seating tolerance. These scooters disassemble in seconds with just four easy steps making transporting from one destination to another a breeze!

BUZZAROUND XL FEATURES:

- 300 lb. weight capacity
- Operating range up to 10 miles on one charge
- 12 Amp batteries standard
- 3.5-inch ground clearance
- LED Angle-adjustable Headlight
- USB Port on the Tiller
- Phone & Cup Holders
- Quick & Easy Disassembly



GB121B-STD

BUZZAROUND SERIES



GB124A-STD

Upgrade with NEW XR Extended Range 22 AH Battery Option!*


The XLSHD offers compact portability with a heavy-duty weight capacity! It features a robust 325 lb. weight capacity and our Comfort Spring Suspension for a smoother ride. It also features a ground clearance of 4" and an operating range of up to 9 miles on the standard 12 AH batteries.

BUZZAROUND XLSHD FEATURES:

- 325 lb. weight capacity
- Comfort spring suspension
- Manually adjustable rear suspension
- Operating range up to 9 miles on one charge
- 12 Amp batteries standard
- 4-inch ground clearance
- LED Angle-adjustable Headlight
- USB Port on the Tiller
- Phone & Cup Holders
- Quick & Easy Disassembly



GB121B-SHZ



GB124A-SHZ

**Packaged separately as an accessory, but ships with unit.*

BUZZ AROUND *Carry On*[®] *Carry On with Life*

**ALL-NEW
Sporty Foldable Scooter!
Airline Friendly!*USB Port Standard!**

Forbes

HEALTH

**NAMES THE BUZZAROUND
CARRYON[®] AMONG
THE BEST LIGHTWEIGHT
MOBILITY SCOOTERS
OF 2022!**



GOLDEN



K. Rippens | Google Review

"I am very happy with my CarryOn scooter. It suits all my priorities: portable, lightweight, comfortable, runs on gravel and grass, easy to fold up for storage. The adjustable tiller is fantastic for bringing the controls within easy reach. It is also a lot of fun to drive."

BUZZAROUND SERIES



The Buzzaround CarryOn foldable scooter is small in stature and big on fun! This super portable, sporty new scooter will take you where you want to go in comfort and style! It features an infinitely adjustable tiller, full front and rear light package, lithium-ion battery (travel battery sold separately*), and 300 lb. weight capacity. The heaviest piece is 47.5 lbs. (without seat and batteries) and the height when folded is only 12". It disassembles easily into two pieces and is designed to store in the smallest compartments. Pack your dreams and go anywhere with the CarryOn!

BUZZAROUND CARRYON FEATURES:

- Designed for minimal storage space
- 4-point carry handles for easy transport
- 12" height when folded (*without seat*)
- 4" ground clearance
- 300 lb. weight capacity
- Flip up armrests
- Infinitely adjustable tiller
- USB port & water bottle holder on tiller
- Disassembles easily
- Long range coverage up to 18 miles
- Easy access 15 AH lithium-ion battery standard
- Optional airline friendly 6.5 AH lithium-ion battery removable for carry on

BUZZ
AROUND
CarryOn[®]

COMPACT TRAVEL



* Optional 6.5 AH lithium-ion battery comes with a stylish bag and can be removed and taken as a carry on while traveling by air, sea or train.

SCAN FOR VIDEO



GB120

LiteRider[®] GOLDEN

LARGE PORTABLE

The unique tubular design of the LiteRider scooters gets you wherever you need to go in a style all their own. The only scooters to come with two cargo baskets standard, the LiteRiders are easy to disassemble, and maneuver smoothly indoors and out.



GL111 & GL141 FEATURES:

- 300 lb. weight capacity
- Wraparound Delta Tiller standard
- Easy to read control panel with backlit battery gauge
- Two storage baskets standard
- Comfortable 17"x16" stadium style seat standard
- Adjustable armrests
- Optional 20" x 18" stadium style seat
- 2-22 AH Batteries; increased operating range up to 12 miles
- Completely wireless disassembly
- 9" front and rear tires
- New Ultrabright LED adjustable headlight and rugged plastic basket standard!

GOLDEN COMPANION

FULL-SIZE LUXURY

Rugged Steel Frame with Weight Capacity up to 400 lbs.!

The Golden Companion scooters are your ticket to ride first class every day! Allow the luxurious, full-size Companions to enhance your life and take you to the places you've always had trouble getting to. Shop in the mall, attend sporting events and enjoy the great outdoors with the unmatched comfort of your favorite new Companion!

ALL COMPANIONS FEATURE:

- All new design for easy disassembly
- Wraparound Delta Tiller standard
- Adjustable LED headlight
- EZ Reach tiller adjustment lever
- Auto lock tiller
- Easy to read control panel
- Easy access charging port
- Black, non-marking, low profile tires on sturdy steel rims
- LED battery gauge
- Comfortable, fully adjustable padded armrests
- Under floor mat troubleshooting guide
- Removable shopping basket with handle
- Four mobility transportation tie downs



Companion GC340 shown with optional power lift seat.



SCAN FOR VIDEO

golden **Avenger**TM

FULL-SIZE OUTDOOR

When you want a rugged, heavy duty scooter to tackle the great outdoors, look no further than the Golden AvengerTM. It features oversized tires with sporty mag wheels, a full-size captain's seat for all day riding comfort, and an impressive operating range of up to 18 miles.



GA541
(Freight Supplement applies.)



AVENGER FEATURES:

- 500 lb. weight capacity
- Adjustable tiller for riding comfort
- Deluxe, reclining high back captain's seat
- Large, 13" front and rear tires
- Adjustable front headlight and four-way flashers
- Large front cargo basket standard
- Operating range up to 18 miles per charge
- Batteries not included: 2-NF22 batteries optional

Golden Patriot™

FULL-SIZE OUTDOOR

If you're looking for the ultimate outdoor driving experience, look no further than the Golden Patriot™. The Golden Patriot's rugged outdoor design features an off-road package, four-wheel independent suspension, a full lighting package and a full-sized captain's seat. With a 400 lb. weight capacity, 13" pneumatic tires, and a battery range of 16 miles, the Golden Patriot will give you the ride of your dreams.



GR575

(Freight Supplement applies.)



PATRIOT FEATURES:

- Full front and rear lighting package with tail lights, headlights and flashers
- Four-wheel independent suspension for smoothly traversing outdoor terrain
- Rugged rubber mat and anti-skid material on the floor deck
- Convenient water bottle holders for two cold beverage bottles
- Full control panel featuring safety horn, flashers, lights, turn signals, speed knob and illuminated battery gauge
- 13" black, pneumatic tires with brushed alloy wheels
- Batteries not included: 2-NF22 or 2-24 batteries optional

Live the Golden Life

Test drive the industry's most comprehensive line of consumer power wheelchairs today!



ONLY
36 lbs.!



GOLDEN
★★★★★

Cheri H. | Google Review

"I got a LiteRider Envy. I have COPD and have issues walking any short distance. This chair enables me to take it with me easily anywhere I need to go. Highly recommend."

The NEW portable LiteRider Envy's design makes our power wheelchairs the easiest to transport!

The LiteRider Envy and NEW Envy LT Power Wheelchairs are both stylish and durable with a tubular design that is unique to the Golden LiteRider family! They feature a large adjustable footplate and an impressive turning radius of 28.5", making them perfect for indoor or outdoor use. They also offer convenient offboard charging and easy disassembly into four pieces. When you are looking for a compact power wheelchair that is fun to drive and easy to transport - look no further than the Golden LiteRider Envy or NEW Envy LT!

The NEW Envy LT features a compact controller and well-styled joystick designed to make everyday tasks easier such as driving up to a table or maneuvering through a doorway. It features an operating range of up to 7.25 miles on two 12V 18AH batteries. The Envy LT offers excellent value and performance for your active lifestyle!



GP161
Available in
Silver Only

LITERIDER ENVY LT FEATURES:

- 250 lb. weight capacity
- 40 Amp joystick controller
- Operating range of up to 7.25 miles
- Adjustable footplate
- Lightweight & Convenient
- Silver Powder Coated Frame
- Offboard charging capability
- Quick & easy disassembly

**Golden's RETAIL
CASH POWER
WHEELCHAIRS!**

SCAN FOR VIDEO



The Envy offers a Dynamic LiNX controller that comes standard and the traditional style makes driving simple and reliable. The Envy, featuring two 22 Amp hour batteries for an operating range of up to 15.5 miles, is Golden's Retail Cash Power Wheelchair!



GP162
Shown with NEW
high back premium
seat option

LITERIDER ENVY FEATURES:

- 300 lb. weight capacity
- Dynamic LiNX controller
- Lightweight & Convenient
- Operating range up to 15.5 miles
- Premium 17" Seat Option
- Disassembles into 4 pieces
- 22 Amp Hour battery for all-day range
- Off board battery charging for on the go

The Envy GP162 is available in 6 Vibrant Colors with a Larger Foot Plate and New Shroud Design!

ENVY GREEN



CHERRY RED



PARADISE BLUE



SUNBURST YELLOW



SATIN SILVER



Experience Constant Contact

This unique suspension system keeps all six wheels on the ground simultaneously, eliminating the rocking sensation when traversing inclines, giving both Compass™ models unparalleled curb and threshold climbing ability. No other power base can match the constant wheel to ground contact of a Compass!

Compass™ models provide superior forward stability making them the safest power wheelchairs to mount or dismount. The Compass line has fewer moving parts than a power base with springs for the best forward stability and no tipping!

golden COMPASS™ SPORT

COMPASS SPORT FEATURES:

- 300 lb. Weight Capacity
- Dynamic LiNX controller
- EZ Lock compatible
- Solid color injection molded plastics
- Nylon drive wheels for reduced weight and modern look
- Silver powder coated arm beams
- Four reflective stickers
- Larger footplate
- Optional 50 Ah batteries for up to 34-mile operating range



Inclines and thresholds are no match for the Golden Compass Sport's true center wheel drive design with our patented springless articulating chassis!

golden COMPASS™ HD

COMPASS HD FEATURES:

- 450 lb. Weight Capacity
- All Black Seat
- Medicare Code: K0825
- Turning Radius: 20"
- Four Post Seat Mount with Slide Seat Bracket
- 75 Amp Controller Standard
- Batteries not included: 2-NF22 batteries optional



The Compass HD has the most stable base in its class!

BUZZABOUT™

Enjoy full-time luxury on the go in the new BuzzAbout! Perfect for travel, the BuzzAbout power wheelchair is in a class all its own in our popular BuzzAround series of mobility products. The BuzzAbout has a stylish captain's seat with contoured support and higher back for longer seating tolerance and comfort. It features motion sensing LED running lights on the footplate and rear handle activated when the joystick is moved. It also offers convenient offboard charging and easy disassembly for transport. The BuzzAbout offers excellent performance and maneuverability indoors and out for your active lifestyle!



**Unique Motion Sensing LED Light
on Footplate & Rear Handle!**

GP164 BUZZABOUT FEATURES:

- 300 lb. weight capacity
- Motion sensing LED running lights on footplate & rear handle
- Quick & easy disassembly into four pieces
- Offboard charging capability
- Dynamic LiNX Controller
- Turning Radius: 28.5"



GP164

TRAVEL LIGHTER

Ride in Style & Secure Your Personal Items with Scooter & Power Wheelchair Accessories from Golden!

TRAVEL CASE FEATURES:

- Stylish, Curved Outer Shell
- Durable Construction
- Black Mesh Partition
- Adjustable Buckle
- Measures 13.77" W x 3.14" H x 9.84" D
- Fits Any Golden Scooter or Power Wheelchair



 GOLDEN



Nancy W. | Google Review

"I love my scooter. It has enabled me to go for "walks" with my family again and enjoy the outdoors. Thanks Golden."

MOBILITY ACCESSORIES

Seat Belts Are Available in 2 Sizes for Added Safety!

60" or 78" Seat Belts fit all Golden Scooter & Power Wheelchair Models



Travel Case



Cane Holder



Scooter & Power Wheelchair Cover



Quad Cane Holder



Rear Basket



Crutch Holder



Walker Holder



Oxygen Tank Holder



Power Elevating Seat*



Optional BuzzAround Shroud Cover Colors:
Pink, Silver & Yellow.
*(Red & Blue standard on most models.
Available for EX and LX models only.)*



Safety Flag

(Included with most models.)



Cell Phone Holder



Cup Holder



Pack'n'Go Armbag

PACK BUNDLE



Rear View Mirror



Y Cable Charging Harness

*Only available on Companion models.

SPECIFICATIONS*

Item	Description	Length	Width	Front Tire	Rear Tire	Center Tire	Casters	Weight Capacity (lbs.) ¹	Maximum Speed (mph)	Ground Clearance	Operating Range ¹ (miles)	Turning Radius	Batteries (2)**	Weight w/ Batteries	Front Half Weight	Rear Half Weight	Seat Weight w/ Arms	Heaviest Piece	Base Weight
GB121B-STD	Buzz XL (3-wheel)	41"	20"	7.5"	8"	N/A	N/A	300	4	3.5"	10	37.5"	12 AH	93	32	24	21	32	N/A
GB121B-SHZ	Buzz XLSHD (3-wheel)	41"	20"	7.5"	8"	N/A	N/A	325	4	4"	9	37.5"	12 AH	100	35	28	21	35	N/A
GB124A-STD	Buzz XL (4-wheel)	40"	20"	7.5"	8"	N/A	N/A	300	4	3.5"	10	45.5"	12 AH	98	37	24	21	37	N/A
GB124A-SHZ	Buzz XLSHD (4-wheel)	40"	20"	7.5"	8"	N/A	N/A	325	4	4"	9	44"	12 AH	105	40	28	21	40	N/A
GB118	Buzz EX (3-wheel)	45"	22"	9"	9"	N/A	N/A	350	5	4"	18	40.5"	U1	154	46	36	23	46	N/A
GB148	Buzz EX (4-wheel)	44"	22"	9"	9"	N/A	N/A	350	5	4"	18	54"	U1	161	53	36	23	53	N/A
GB119	Buzz LX (3-wheel)	45"	22"	9"	9"	N/A	N/A	375	5	4"	18	40.5"	U1	168	50	36	32	50	N/A
GB149	Buzz LX (4-wheel)	44"	22"	9"	9"	N/A	N/A	375	5	4"	18	54"	U1	174	57	36	32	57	N/A
GB120	BuzzCarryOn (4-wheel)	38"	19.5"	6"	7"	N/A	N/A	300	3.6	4"	18	48.75"	15 AH	65.5	47.5	N/A	12	47.5	N/A
GL111	LiteRider (3-wheel)	41.5"	22"	9"	9"	N/A	N/A	300	5	3.75"	12	39"	22 AH	125	44	28	21	44	N/A
GL141	LiteRider (4-wheel)	43.25"	22"	9"	9"	N/A	N/A	300	5	3.75"	12.5	48"	22 AH	133	52	28	21	52	N/A
GC240	Companion (3-wheel)	42.25"	23.5"	9"	10.75"	N/A	N/A	350	5.3	3.75"	19.9	37"	U1	163	45	48	26	48	N/A
GC340	Companion (3-wheel)	45.5"	23.5"	9"	10.75"	N/A	N/A	400	5.3	4.5"	20.9	41"	U1	176	56	48	26	56	N/A
GC440	Companion (4-wheel)	47.5"	23.5"	10.75"	10.75"	N/A	N/A	400	5.3	4.5"	18.9	56.5"	U1	194	74	48	26	74	N/A
GA541	Avenger	57.5"	24.5"	13"	13"	N/A	N/A	500	7	4"	18 ²	67"	NF22**	265	68	68	55	68	N/A
GR575	Patriot	57"	28"	13"	13"	N/A	N/A	400	10	6.75"	16	76.5"	75 AH**	302	N/A	N/A	44	152	152
GP161	LiteRider Envy LT	33.25"	22.5"	N/A	8"	N/A	6"	250	3.5	2.5"	7.25	28.5"	18 AH	113	34.5	30.5	23	34.5	N/A
GP162	LiteRider Envy	36"	22.5"	N/A	9"	N/A	6"	300	3.5	2.5"	15.5	28.5"	22 AH	123	37	36	23	37	N/A
GP164	BuzzAbout	36"	23"	N/A	9"	N/A	6"	300	4	2.5"	11	28.5"	22 AH	124	29	33	34	34	N/A
GP605	Compass Sport	39.5"	24"	N/A	N/A	10"	6"	300	4	3.25"	24	21"	U1	179	N/A	N/A	40	97	97
GP620	Compass HD	40"	25"	N/A	N/A	10"	6"	450	4	3.5"	12 ²	20"	NF22**	226	N/A	N/A	51	101	101

1. Battery range at 200 lbs., and will vary due to rider weight, drive surface, terrain and battery type.

2. Battery range at maximum weight capacity, and will vary due to rider weight, drive surface, terrain and battery type.

Brakes: All Golden Technologies scooters come standard with Dynamic/Electromechanical braking systems.

Rear wheel drive: All Golden Technologies scooters come standard with a sealed transaxle, 24 volt motor.

* All weight shown in lbs. ** Batteries not included.

All products feature Freewheel Mode, Electronic Speed Control and Off-Board Charger.



<https://m.youtube.com/user/GoldenTechChannel>



Complies with
UL 73 and
CSA C22.2 No. 68



www.goldentech.com



www.facebook.com/goldentechnologies

Note: Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

Visit your Authorized Golden Dealer: